ZENET

MASSAGER ZET-719

MASSAGER ZET-719
USER MANUAL
Voltage / frequency 220V / 50-60Hz
Power 25W

Dear customer,

Thank you for buying one of
Models of ZENET massagers.

Please read the operating instructions and keep them
for future reference.

PRINCIPLE OF THE MASSAGE:

Massage is one of the treatments in traditional Chinese medicine. With the help of the massage one can act on certain parts of the body to treat a certain illness.

Massage is an excellent therapy, with whose help you can make your body stronger and stronger. Massage is widely used in modern medicine. According to the theory of Chinese medicine, the surface of the human body is directly connected to the internal organs with the help of certain points on the body. If something is wrong with the internal organs, the disease affects the surface of the body.

The massage of the main and auxiliary points on the body surface also stimulates the blood circulation.

Massage promotes the body's resistance, the absorption of the inflammatory process and can relieve muscle cramps and pain.

DESCRIPTION OF THE DEVICE:

1. The power regulator
2. The lever
3. The straps
4. The massage nozzle
5. The delimiter
6. The function selection

APPLICATION:

Make sure that the device is in the “OFF” state and plug the power supply unit

unit into the corresponding voltage socket. Set the switch to the «M» position, the device is switched on.

If you turn the speed control clockwise, the massage intensity gradually increases. The reverse action leads to a reduction in intensity.

The duration of the massage should not exceed 15 minutes in order not to cause excessive stimulation. At the end of the massage, switch off the power supply and pull out the plug.

GUARANTEE LIABILITIES:

1. Guarantee of the device according to the operating rules - 1 year from the date of sale.
2. During the warranty period, the device must be replaced if any defects caused by the manufacturer are found. The sales date and the stamp of the trade organization are required.
3. The guarantee does not apply to devices with:
 - mechanical damage;
 - defects resulting from a violation of the terms of use or incorrect actions by the owner of the device;
 - deficiencies due to natural disasters;
 - defects caused by the ingress of foreign bodies and liquids as a result of self-repair or self-introduction of structural changes;
 - used for commercial or industrial purposes.

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