ZENET

MASSAGER ZET-709

MASSAGER ZET-709
INSTRUCTION MANUAL
Voltage / frequency 220V / 50-60HZ
Power 25W

Dear customer,
Thank you for buying one of
Models of ZENET massagers.

Please read the operating instructions and keep them
for future reference.

DESCRIPTION:

An advanced manual massager that provides vibration massage for your back and body, for both acupressure and fingers!
The ergonomic shape of this massager makes it possible to massage the entire body with little effort: collar area, shoulders, back, legs and feet.
Safe interchangeable silicone nozzles (3 types of 2 pieces). Here you can select the type of massage and its intensity.
Nozzles perform various functions: superficial and deeply penetrating muscle massage to prevent osteochondrosis or to treat cellulite massage to accelerate subcutaneous fat. The presence of a thermal effect allows you to treat sinusitis, runny nose, hematoma and varicose veins.
The positive effect of this hand massager has been proven in practice for many years.

FUNCTIONS:

- Elimination muscle fatigue and rheumatism, shoulder cramps.
- Restoration of blood circulation and partial fat burning.
- Thermal exposure to improve exposure.

OPERATING MANUAL:

1. Insert the plug into an electrical outlet.
2. Press the «ON / OFF» key to start or stop the instrument.
3. Press the «SPEED» key to increase or decrease the strength of the massage (4 speeds).
4. Press the «PROGRAMS» key to select one of three massage modes.

GUARANTEE LIABILITIES:

1. Guarantee of the device according to the operating rules - 1 year from the date of sale.
2. During the warranty period, the device must be replaced if any defects caused by the manufacturer are found. The sales date and the stamp of the trade organization are required.
3. The guarantee does not apply to devices with:
 - mechanical damage;
 - defects resulting from a violation of the terms of use or incorrect actions by the owner of the device;
 - deficiencies due to natural disasters;
 - defects caused by the ingress of foreign bodies and liquids as a result of self-repair or self-introduction of structural changes;
 - used for commercial or industrial purposes.

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